



COMPUTER SPECIFICATIONS

ITEM	Recommended
Processors	<ul style="list-style-type: none">• Intel i5 or better (NO Small Form Factor (SFF) or Micro computers) *• NO VIA or SIS chipsets• Regular or Active USB extension cable to connect sensors to the pc. <p><i>*SFF pcs are usually low powered and lack space to add additional cards if needed.</i></p>
Operating System	<p>Workstation:</p> <ul style="list-style-type: none">• Windows 10, Windows 11 – Professional versions* <p>Windows Server Versions:</p> <ul style="list-style-type: none">• 2012 / 2016 / 2019 <p>Check compatibility of hardware drivers with version of Windows.</p> <p><i>*Windows 7 / Server 2008 Microsoft Support ended January 2020</i></p>
Random Access Memory (RAM)	<ul style="list-style-type: none">• Workstation: 8 GB or better• Server: 16 GB or better <p>More RAM increases performance</p>
Hard Drives	<ul style="list-style-type: none">• Workstation: 80+ GB• Server: 500+ GB for data storage
Monitor and Video Card	<ul style="list-style-type: none">• Widescreen monitor (1366 x 768 minimum resolution)• Dedicated Video Card with 1GB+ Memory
Network and Support	<ul style="list-style-type: none">• 1 Gigabit WIRED Network Card (Wireless NOT recommended)• HIGH SPEED Internet connection for support
Backup System	<ul style="list-style-type: none">• Windows approved backup system for CADI Patients' folders.• Backups should be stored off site.